# **ENDO STRATUS®**

Irrigation Pump and CO<sub>2</sub> Insufflator | **Procedure** 







#### THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, only Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

# Advanced CO<sub>2</sub> Insufflation for Faster Procedures and Better Patient Outcomes.

## Quicker Patient Turnaround\*

- Improved cecal intubation rates<sup>1</sup>
- · Greater small bowel intubation depths2
- 38% reduction in nursing attention

## Patient Pain and Discomfort CO2 vs Room Air Co2 Room Air 6 5 1 During 1 Hour After 3 Hours 6 Hours

## **Higher Patient Satisfaction**

- · Less discomfort, pain and bloating<sup>3,4,1</sup>
- · Quicker recovery times4
- · Decreased gas distension
- Absorbed 150 times faster than room air. and is promptly eliminated via the lungs

Warm water technology minimizes chances of spastic colon Heated CO<sub>2</sub> – 98.6 degrees

lens cleaning and insufflation

Standard male luer output -Compatible with ENDO SMARTCAP™

Easy push button ON/OFF operation

Variable flow rates – low, medium & high

- 1. Singh, R., Neo, E., & Nordeen, N. (2012, July). Carbon dioxide insufflation during colonoscopy in deeply sedated patients. World Journal of Gastroenterology, 18(25), 3250-3253. 2. Bretthauer, MD, PhD, M. (2010). Turning science into clinical practice - the case of carbon dioxide insufflation. Endoscopy 2010, 42(), 1104-1105. doi:http://dx.doi.org/10.1055/s-0030-1255973
- 3. Bretthauer, MD, PhD, M. (2007, September). Carbon dioxide insufflation improves intubation depth in double—balloon enteroscopy: a randomized, controlled, double—blind trial. Endoscopy 2007, 39(), 1064-1067. doi:DOI 0.1055/s—2007—966990 4. Lynch, MBA, BSN, RN, Hayes BSN, RN, CGRN, Buffurn DNSc, APRN, CS. (2012). Insufflation Using CO2 vs Room Air During Colonoscopy: Comparison of Patient Comfort, Recovery Time, and Nursing Resources [PowerPoint slides]. Veterans Affairs Medical Center San Francisco.



#### **PROCEDURE**

Reducing the risk of patient cross contamination is at the forefront of infection prevention. Cantel innovates infection control products designed to improve patient outcomes, while increasing procedural efficiency.

# **ENDO STRATUS®**

Irrigation Pump and CO<sub>2</sub> Insufflator | Procedure



## ENDO STRATUS® CO<sub>2</sub> Insufflator

Compatible with

- Tank (D or E) or DISS Wall CO<sub>2</sub>
- · All major GI Endoscopes

## ENDO STRATUS® Irrigation Pump

Adjustable water bottle holder

Compatible with

· All major GI Endoscopes

Pump includes comfortable, universal foot pedal

Warm Water Irrigation Improves
Your Visibility and Your Patient's Comfort

Easy and quick to adjust flow rate

Compatible with MEDIVATORS ENDOGATOR® tubing and connectors

Automatic prime button provides instant irrigation upon foot pedal depression



# 50.05" (127.1 cm) 32.57 (82.7)16" (40.6 CM)

## CO<sub>2</sub> Carts

- · Conveniently transport and store two CO<sub>2</sub> tanks.
- ENDO STRATUS<sup>®</sup> Irrigation Pump and CO<sub>2</sub> Insufflation units stack neatly on top platform.
- Non-porous straps secure ENDO STRATUS Irrigation Pump and CO<sub>2</sub> Insufflation units for easy transportation between rooms.
- · Basket provides additional storage.

ENDO STRATUS® Carts Ordering Information

Part Number	Description
EGA-7030	ENDO STRATUS IV Style Cart
EGA-7032	ENDO STRATUS IV Style Cart with Dual Yoke
EGA-7031	ENDO STRATUS Box Style Cart
EGA-7033	ENDO STRATUS Box Style Cart with Dual Yoke



17.80"(45.2 CM)



### ENDO STRATUS® Pumps Ordering Information

Description
ENDO STRATUS Irrigation Pump
ENDO STRATUS Insufflator Unit
ENDO STRATUS Irrigation Pump & CO <sub>2</sub> Insufflator
Tank Yoke
Dual Tank Yoke
3 Foot Wall Hose (Low Pressure)
9 Foot Wall Hose (Low Pressure)
3 Foot Tank Hose (High Pressure)
9 Foot Tank Hose (High Pressure)

ENDO STRATUS®, ENDO SMARTCAP® and ENDOGATOR® are registered trademarks of Medivators Inc.





**Global Headquarters** 14605 28th Avenue North Minneapolis, MN 55447-4822 Phone: +1.763.553.3300 Fax: +1.763.553.3387

Call: 1.800.328.3340 (Customer Service) **TO PLACE AN ORDER** 

Fax: 1.800.686.8493

Email: custserv@medivators.com